



P R O M O T I O N

O C T - D E C 2 0 2 3

PRIME LINE



AeroAcoustic blade



Glass rod fiber optics



Reverse flow protection



Automatic pressure regulation¹



Refined ceramic bearings



Scratch resistance
Iconic Prime body

**25
W**

**Powerful & Precise
Power Head - 25W**

**55
dB**

**Incredibly Quiet
Only 55 dB noise**

Available for
KaVo*, NSK*, W&H*
and Sirona* Couplings



**BUY 2
&
GET 1
FREE**

Coupling Connection

Kavo
NSK
W&H
Sirona

Power Head 25 Watt

HP21KL
HP21NL
HP21WL
HP21SL

Small Head 18 Watt

HP22KL
HP22NL
HP22WL
- N/A -

\$ 1,683.00



**Now with 2 Years
FairCare Warranty²**

Unconditional warranty
except physical damage
on Prime Line High-speeds

COUPLINGS



✓ MK-dent & Kavo
with LED & water
• **QC6016KW**
\$ 495.00



✓ NSK connection
with LED
• **QC6016NT**
\$ 385.00



✓ W&H connection
with fibre optic
• **QC5016W**
\$ 385.00



✓ Sirona connection
with LED & water
• **QC6016SW**
\$ 495.00

✓ Non-optic
• **QC4014K**
\$ 296.00

✓ Non-optic
• **QC4014NT**
\$ 296.00

✓ Non-optic
• **QC4014W**
\$ 296.00

✓ Non-optic
-NA-

1) Only available for Kavo coupling 2) FairCare warranty applicable on Prime Line high speeds sold after 1st Oct 2018

PRIME LINE CONTRA - ANGLES ♡

- Buy 3 - Get 15% OFF
 - Buy 4 - Get 20% OFF
- Mix & Match



2 Years FairCare Warranty
Unconditional warranty
except physical damage



⚙️ 1:5 Red ●

⌚ Max. 200,000 rpm

○ Internal water

🔧 For FG-burs
(Ø 1.6 mm)

● LP15L (F/Optic)

RRP: \$ 1,826.00

Buy 3 - \$1,552.10 ea

Buy 4 - \$1,460.80 ea

⚙️ 1:1 Blue ●

⌚ Max. 40,000 rpm

○ Internal water

🔧 For CA-burs
(Ø 2.35 mm)

● LP11L (F/Optic)

RRP: \$ 1,287.00

Buy 3 - \$1,093.95 ea

Buy 4 - \$1,029.60 ea

⚙️ 5.4:1 Green ●

⌚ Max. 7,407 rpm

○ Internal water

🔧 For CA-burs
(Ø 2.35 mm)

● LP54L (F/Optic)

RRP: \$ 1,320.00

Buy 3 - \$1,122.00 ea

Buy 4 - \$1,056.00 ea

⚙️ 1:1 Blue ●

⌚ Max. 40,000 rpm

○ Internal water

🔧 For HP-burs
(Ø 2.35 mm)

● LP01L (F/Optic)

RRP: \$ 1,045.00

Buy 3 - \$888.25 ea

Buy 4 - \$836.00 ea

AIR MOTORS ♡



- ✓ Internal water
- ✓ LED

● AM1006
\$ 990.00



- ✓ Internal water

● AM1004
\$ 748.00



- ✓ External water

● AM0014
\$ 440.00



● AM0004
\$ 440.00

NEW

KAVO COMPATIBLE HEADS ♡



- Kavo 68LH compatible
- Blue (1:1)
- LCH11
\$ 440.00



- Kavo 67LH compatible
- Green (2:1)
- LCH21
\$ 440.00



- Kavo 31LR compatible
- Green (2:1)
- LCH21P
\$ 440.00

ECO LINE TURBINES



Now with 1 Year FairCare Warranty¹
Unconditional warranty except physical damage



Now with 24W*

- 62 dB Operating Noise
- AeroPerformance blade
- Refined ceramic bearings
- Glass rod fiber optics
- Smooth Eco body



Available for
KaVo*, NSK*, W&H*
and Sirona* Couplings



**BUY 3
or MORE**
Only
\$759
each

**3 x Eco Line HS
with F/O for any
coupling connection**

**Savings
from \$392**

Coupling Connection	Power Head 24W (F/Optic)	Small Head 16W (F/Optic)	Power Head 24W (N/Optic)	Small Head 16W (N/Optic)
KaVo	HE21KL	HE22KL	HE21K	HE22K
NSK	HE21NL	HE22NL	HE21N	HE22N
W&H	HE21WL	HE22WL	HE21W	HE22W
Sirona	HE21SL	- N/A -	- N/A -	- N/A -

• \$ 890.00 (F/Optic)

• \$ 649.00 (N/Optic)

Buy 3 - Get 10% OFF

ECO LINE POWER TURBINE

NEW - With 28w POWER



• HE20KL
\$1,265.00

Buy 3 or more
Only **\$1,089 ea**

18
MONTHS
WARRANTY

Available for only KaVo coupling

60 dB Operating Noise

AeroPerformance blade

Refined ceramic bearings

Glass rod fiber optics

Smooth Eco body

ECO LINE FIXED CONNECTION TURBINES

- Direct 4 hole connection
- Refine ceramic bearings
- Push Button System
- HE17 - Standard head with 4 hole water spray
- HE18 - Small head with 3 hole water spray



\$ 539.00

Buy 3 - Get 10% OFF




1) FairCare warranty (except physical damage) applicable on Eco Line sold after 1st Jan 2019 2) Eco Line Power Head 24 Watt

• Buy 3 - Get 15% OFF
 • Buy 4 - Get 20% OFF
Mix & Match



FairCare Warranty¹
 1 year unconditional warranty
 except physical damage






-  1:5 Red ●
-  Internal water
-  For FG-burs
(Ø 1.6 mm)

• LE15L (F/Optic)
 RRP: \$ 1,485.00
Buy 3 - \$1,262.25 ea
Buy 4 - \$1,188.00 ea

• LE15 (Non-optic)
 RRP: \$ 1,287.00
Buy 3 - \$1,093.95 ea
Buy 4 - \$1,029.60 ea






-  1:1 Blue ●
-  Internal water
-  For CA-burs
(Ø 2.35 mm)

• LE11L (F/Optic)
 RRP: \$990.00
Buy 3 - \$ 841.50 ea
Buy 4 - \$ 792.00 ea

• LE11 (Non-optic)
 RRP: \$770.00
Buy 3 - \$ 654.50 ea
Buy 4 - \$ 616.00 ea






-  5.4:1 Green ●
-  Internal water
-  Screw-in & Snap-on
Prophy Cups

• LE54LP (F/Optic)
 RRP: \$ 1,210.00
Buy 3 - \$ 1,028.50 ea
Buy 4 - \$ 968.00 ea

• LE54P (Non-optic)
 RRP: \$935.00
Buy 3 - \$ 794.75 ea
Buy 4 - \$ 748.00 ea





-  1:1 Blue ●
-  Internal water
-  For HP-burs
(Ø 2.35 mm)

• LE01L (F/Optic)
 RRP: \$935.00
Buy 3 - \$ 794.75 ea
Buy 4 - \$ 748.00 ea

• LE01 (Non-optic)
 RRP: \$748.00
Buy 3 - \$ 635.80 ea
Buy 4 - \$ 598.40 ea



**Buy More
 Save More**



-  1:1 Blue ●
-  Push button
- LB11/11
 Buy 1 - \$418.00



Buy 3 - \$339.00 ea
Buy 4 - \$319.00 ea



-  1:1 Blue ●
-  Latch Type
- LB11/11C
 Buy 1 - \$418.00


Buy 3 - \$339.00 ea
Buy 4 - \$319.00 ea



-  4:1 Green ●
-  Snap-on Cups
- LB41/11P
 Buy 1 - \$ 462.00

Buy 3 - \$396.00 ea
Buy 4 - \$385.00 ea



-  1:1 Blue ●
- LB02
 Buy 1 - \$ 418.00

Buy 3 - \$339.00 ea
Buy 4 - \$319.00 ea

¹ FairCare warranty (except physical damage) applicable on Eco Line sold after 1st Jan 2019

\$ Buy 3 - Get 1 FREE

1 Year Warranty



EMS Non-optic Type EN-041

- SC21E
Buy 1 - \$440.00
Buy 3 - Get 1 FREE



EMS Piezon Non-optic Type EN-061

- SC21E2
Buy 1 - \$539.00
Buy 3 - Get 1 FREE



EMS Piezon LED Type EN-060A

- SC21EL
Buy 1 - \$715.00
Buy 3 - Get 1 FREE



KaVo Piezo LED Type 1.007.3995

- SC21KL* **NEW**
Buy 1 - \$715.00
Buy 3 - Get 1 FREE



Acteon-Satelec Non-optic Type F 12 281

- SC21S
Buy 1 - \$ 440.00
Buy 3 - Get 1 FREE



Acteon-Satelec LED Type F 12 609

- SC21SL
Buy 1 - \$715.00
Buy 3 - Get 1 FREE

*Only for use with EMS compatible scaler tips



Buy 3 - Get 1 FREE while stocks last

- ✓ High-quality Diamond-Like-Carbon coating
- ✓ Extremely hard and wear-resistant surface

• **SCALING DLC COATED** Only available for Satelec fitting

ST02SD SATELEC 1S
\$88.00

ST03SD EMS A
\$88.00

ST04SD SATELEC 2
\$132.00

• **PERIO DLC COATED** Only available for Satelec fitting

ST06SD SATELEC 10Z
\$132.00

ST10SD SATELEC 10P
\$154.00

• **ENDO DLC**

ST08SD SATELEC S07
\$154.00

• **SCALING STAINLESS STEEL** Only available for Satelec fitting

ST02S SATELEC 1S
\$63.80

ST03S EMS A
\$63.80

ST04S SATELEC 2
\$88.00

• **PERIO STAINLESS STEEL** Available for EMS & Satelec fittings

ST07E/ST07S EMS DPL3 SATELEC 10X
\$110.00

ST09E/ST09S EMS HPL
\$110.00

ST10E/ST10S SATELEC 10P
\$110.00

ST11E/ST11S EMS RT1
\$110.00

Universal Torque Wrench

TORQ CONTROL®

Precision instrument with a Torque Adjustment System for tightening of Prosthetic components at specific Torque, for any implant system.



Automatic Declutching when Selected Torque is reached

- Compliance with Protocols for screw tightening
- Limits the risk of the screw fracturing or loosening




Universal connection of Prosthetic mandrels

- 7 different Torque levels (from 10 to 35 N.cm)
- Can be used with ALL Implant Systems



Micro-Head and Ergonomic Handle

- Easy access in restricted area with Excellent Visibility



Easy and Quick to use

- Reduced motion compared with using a ratchet wrench while tightening



UNIQUE

\$1,210.00


Lightweight
ONLY 135 gms

Order code: 15501

3 Easy Steps



1. Please select the recommended Torque



2. Start rotating



3. Screw automatically stops when selected torque is reached



www.anthogyr.com

Anthogyr
A Straumann Group Brand

▼ PROPHY LINE AIR POLISHING DEVICE



**Buy 1 Polisher and receive
• 1 bottle of Propphy Powder FREE
and
• 50% discount on 1 additional tip**

• PR1011
\$1,375.00



• Ergonomically shaped, compact and light-weight air-driven polishing device

• 360° rotatable outer handle due to unique RotaGrip® for simple and fatigue-proof operation



• Propphy Powder (300 gms)
\$ 60.50 each

• Lemon - PP1012

Additional Propphy tip



• SP3011
\$ 374.00

Coupling connection order codes:

Kavo - PR1011K, NSK - PR1011N, W&H - PR1011W, Sirona - PR1011S

▼ LEDs



Buy 4 - Get 10% OFF

- ✓ Quality illumination, No Glare
- ✓ 5500K Natural Daylight Colour
- ✓ Clear Vision & Sharper Contrast
- ✓ Longer Durability & Economical



\$88.00

- BU8012 - KaVo & MK-dent coupling
- BU8012KM - KaVo & MK-dent motors
- BU8012WH - W&H handpieces
- BU8012WLT - W&H airmotors LT25
- BU8012WEA - W&H motors 40LT, 50 LT & A51LT
- BU8012SG - Sirona coupling and motors
- BU8012N - NSK coupling
- BU8012NT - NSK motors TiM40 & Ti205L
- BU8012NM - NSK coupling KCL,SCL,PTL LED

▼ PREMIUM LUBRICATION OIL



- ✓ Fully Synthetic
- ✓ For any brand of handpieces and motors

• LU1011 (500 ml)
\$ 68.20



- ✓ Fully Synthetic
- ✓ For KaVo Quattrocare 2 lubrication system

• LU1022 (500 ml)
\$ 68.20



- ✓ For W&H Assistina lubrication system

• LU1031 (500 ml)
\$ 88.00

MANUFACTURER



MK-dent GmbH
Marie-Curie-Str. 2
22941 Bargteheide
GERMANY
www.mk-dent.com

YOUR SPECIALIST DEALER



Suite 201/ 15-17 Park Street,
South Melbourne, VIC 3205

T: 1 800 898 749

E: info@innodent.com.au

W: www.innodent.com.au

• All catalogue prices are inclusive of GST • Promotion valid until 31st Dec 2023 or while stocks last

KaVo, W&H, Sirona, NSK, Bien-Air, EMS and Acteon-Satelec are registered trademarks. MK-dent GmbH has no economic connection with the aforementioned companies.